



WF-20(W)

\$279

Product Images



Description

- Compatibility Details: [WF-20 Compatibility Chart](#)
- How to use WF-20 with DinoCapture: [Support Article](#)

Coupling with Dino-Lite AF series, the WF-20(W) is a Wireless Streamer adding further versatility, flexibility, and agility to Dino-Lite microscope. Thanks to stable WiFi image transmission and long battery life, the WF-20(W) is well suited but not limited to field applications. The WF-20(W) can also be configured as a wireless router for staying connected with the internet or a network during usage. The WF-20(W) Connect allows AF models to be used wirelessly in connection with an iOS (iPhone/iPad) or Android app on any tablet, smartphone, or even in the DinoCapture software for Windows computers. The WiFi streamer is ideal for field work or presentations and can stream from one Dino-Lite to multiple devices. When using the iOS app or DinoCapture, users can perform measurements and take advantage of the automatic magnification reading (AMR) with compatible Dino-Lite AMR and measurement models. The WF-20(W) can be used anywhere to create a WiFi signal for streaming the Dino-Lite live images. The free DinoConnect App for WF-20(W) can be downloaded from the App Store or Google Play.

*For USB models see the [WF-10](#)

Additional Information

Warranty Period	1 year
Manufactured in	Taiwan
Service & Support	U.S. (www.dinolite.us) English & Spanish
Regulatory Approval	CE, FCC
Dimensions	58mm (L) x 36.5mm (W) x 32mm (H)
Weight	57g
Battery Life	2.5 hours (continuous use)
Compatible with	"AF" models only, App: iOS 8.0 or newer Android 4.4 or newer
Package Includes	WF-20(W) Wireless Adapter

