



MSAA105

\$9

Product Images



Description

Curved Lens cap adapter for AD413 series

Additional Information

Warranty Period	1 year
Manufactured in	China
Service & Support	U.S. (www.dinolite.us) English & Spanish
Regulatory Approval	ROHS
Dimensions	30mm(D) x 21mm(L)
Weight	3g
Material	Plastic
Compatible with	"AD" models excluding "Z" & "L"
Package Includes	Lens Cap Adapter
Package Dimensions	7.6cm(L) x 5.1cm(W) x 3.2(H)
Package Weight	27g

