



Why choose Dino-Lite microscopes?

\$0

Product Images



Professional yet affordable

- Developed for professionals yet economical for consumer or casual use
- Award winning product with U.S. design patent
- #1 Selling Handheld Microscope in the World
- If the wide selection of products do not meet your needs, custom options are available

Dino-Lite® Digital Microscope



Answer

Special Features

- Exclusive features such as **Enhanced Depth of Field (EDOF)*†** and **Extended Dynamic Range (EDR)*†**
- **One-Touch Auto Focus** available with EDOF models*
- **Wide range of special lighting models†**
- **Automatic Magnification Reading (AMR)†** for easier and more accurate measurements
- **Axial lighting** models providing darkfield and brightfield viewing abilities*
- **Interchangeable front caps**, provide solutions for handheld viewing distances and lighting control for imaging clarity†
- **Long and Far working distance** models that provide additional focal distance to meet application needs

Hardware

- **Unique design** with single glass lens and variable optical magnification
- **Built-in Polarization** for reducing or eliminating glare†
- **Innovative Flexible LED Control (FLC)** with adjustable intensity and sectional lighting†
- **New modular design** available with option for wireless microscopy via WF-20 adapter to iOS and Android, smartphones, tablets, and WiFi computers
- **DinoDirect** connection to Android smartphones and tablets; **WF-10 WiFi** wireless connection to Apple iPhone and iPad
- **Dino-Lite** series microscopes are the original patent design for handheld microscopes, beware of possible design copies
- **Huge selection of accessories** designed to meet needs of users in many applications
- **Scroll Lock** to maintain calibrated installations†
- **MicroTouch** capacitive touch sensor technology for easy handheld use†

Software

- **Actively maintained, user friendly, advanced DinoCapture 2.0*** software included with free updates and frequently added new features with **release notes** and **online user guide**.
- **Highly accurate measurement tools and calibration features†**
- **Snap-to-Edge†** capability to automatically detect edges when performing measurements
- **Customizable shortcut keys*** to help improve workflow
- **Edge detection control†**, **grid lines†**, **reticle†**, **custom overlay**, **transparency**, and many more useful features
- **Software** can be installed on multiple computers at no additional cost
- **Windows** and **Mac** software
- **Software Development Kit (SDK)** available for free with signed license agreement providing users the ability to customize or integrate the Dino-Lite into a system
- **Customized automated image analysis** through compatibility with **OpenCV** and **Image-Pro Premier**

Service

- **2 Year warranty** on Dino-Lite microscopes and Dino-Eye eyepiece cameras
- **1 Year warranty** on all Accessories designed for Dino-Lite
- **Free US based technical support** providing help with software and hardware questions over the phone, via email, or even remote assistance
- **Complimentary 30 minute online training session** for Dino-Lite and DinoCapture
- **70+ models** to choose from and tailored to meet your application needs
- **Wide selection of industry specific accessories**
- **99% of all items** are in stock and readily available to ship same business day
- ✓ **Available on certain models only.** Please check individual product pages to confirm compatibility.
- * Available on Windows only.
- † Available with specific Edge Series models only.

