





# MS62P

\$99

## Product Images





## Description

---

Rolling rack with MSA2Z holster great for flat surfaces that require viewing in all direction movement. Move the stand easily with 4 rolling ball, clear grip to ensure steady movement. The MS62P can be used in normal settings also as a stand for your Dino-Lite handheld microscope making it another versatile accessory to help you maximize your performance.

## Additional Information

---

Warranty Period	1 year
Manufactured in	China
Service & Support	U.S. (www.dinolite.us) English & Spanish
Regulatory Approval	ROHS
Dimensions	10.1cm (W) x 10.1cm (L) x 4.4cm (H)
Weight	227g
Material	Stand: Metal, Holster: Plastic
Flexibility	Multi Angle
Compatible with	All Dino-Lite Handheld Microscopes
Package Includes	Stand, Holster
Package Dimensions	15.6cm(L) x12.1cm(W) x8.6cm(H)
Package Weight	213g

